

# What are my options?



**Digital** 

**Paper** 

**Print edition** 

(Paper Manuals and all Record Forms)




If you use a desktop or laptop computer or have a tablet (not an iPad), do you have web access?



If you have access to iPad®s and you have web access

No

**Option 1** (Requires at least 1 iPad)





Q-global°

Web-based

reporting



Q-interactive®

# administration, scoring, &

Yes

**Purchase** 

WRAT5 Digital Kit on Q-global, which includes:

Manual

Spelling Card

Norms Book

- Reading Card
- Sentence Comprehension Card Set

• Pay-per-use administration, scoring, and reporting combined

**Digital** 

administration, web-based

scoring, & reporting

## **Administer**

Print needed pages or download the Administration Directions document from the Q-global Resource Library and open the Digital Stimulus Book to display the visual stimuli on your screen (an iPad or any single device greater than 9.7" diagonal).

Record examinee responses on a paper Record Form.

#### Score

Score online, entering either item scores or raw scores from the paper Record Forms (scoring purchased separately). Display pages from your digital Manual or Technical Manual on your screen, and use search, highlight, or bookmark functions to find information you use frequently.

#### **Administer**

**Purchase** 

Annual license

Use two iPads-one for the examinee to see the test stimuli and one for you to capture the examinee's responses. Responses to the receptive tests are scored automatically.

#### Score

Automated scoring included. View scaled scores when you complete testing. View index scores after each test session on Q-Central.

or scoring options are available without web access at this time.

No digital test

administration

### Pearsonclinical.co.uk

Copyright © 2019 Pearson Education, Inc. or its affiliates. All rights reserved. Pearson, Q-interactive, and Q-global are trademarks, in the U.S. and/or other countries, of Pearson plc. Wide Range and WRAT are trademarks of Wide Range, Inc. CLINA15775-23040 ML 2/20

